

# **Updates for Raith Systems**

Raith's Electron Beam Lithography and Focused Ion Beam systems are reliable and long lasting. These systems are an efficient and essential part of your tasks and processes.

Raith offers a huge variety of products to keep your Raith system up to date and reliable for many years. So your Raith system will remain the ideal solution for the jobs that lie ahead!



# Not too late yet to get Win10 (64bit)

Microsoft has stopped support of Windows XP and Windows 7. Bug fixing and crucial security updates are no longer available. If you have not already done so, Raith strongly recommends you to update your system to a Windows 10 environment. Raith supplies preconfigured update packages to avoid any delays during the installation and keep downtime to a minimum. Please contact your Raith Service office to get your quote.

## Good reasons, no more risks:

- · Windows 10 security updates will keep your system safer
- Avoid local IT and network restrictions and ensure Raith Remote Support remains available
- · Required for updates of Raith Nano Suite

## Speed up - new aperture driver

A new type of aperture driver is now available for many Raith systems. Benefit from smoother and faster usage by updating the aperture driver of your aging system. The new aperture driver has a step size 2 times better than the old version (1 $\mu$ m vs. 2  $\mu$ m). Causing improved gun alignments you will benefit from a better beam quality. Additionally, changing between user modes will speed up.

# The new aperture driver is faster with a higher precision!

Contact your Raith service center to check its availability for your Raith system!

Contact us! support@raith.com Phone: +49 231 95004 499

### Service contact

#### Europe

#### Germany / Head office

Raith GmbH Konrad-Adenauer-Allee 8 44263 Dortmund, Germany

Phone +49 231 95004 499 Email support@raith.com

#### Netherlands

Raith B.V.
De Dintel 27a
5684 PS Best, The Netherlands
Phone +31 499 336 888
tsc.helpdesk@raith.com

#### **Americas**

Raith America Inc., Islandia, NY
Phone +1 855 558 9558
Email support@raithamerica.com

#### Asia / Pacific

Raith Asia Ltd., Hong Kong Phone +852 2887 6828 Email support@raithasia.com

#### China

Raith China Co., Ltd., Beijing
Phone +86 10 828679 22
Email support@raithchina.com

### India

Raith India Pvt. Ltd., Bangalore Phone +91 80 2838 4949 Email support@raithindia.com



For further contact information, please visit our website:

www.raith.com

06-2020

